


# The iFix™ System

## PEEK™ INTERFERENCE SCREW



	iFix™	Metal Screws	Bio Screws
Strong, Rigid Fixation	✓	✓	
Radiolucent	✓		✓
Bio-Inert	✓	✓	
Easily Revised	✓		✓
No Tapping Required	✓	✓	
Ideal Modulus	✓		

**FAST. STRONG. COMPLETE.**



# OPTIMAL PERFORMANCE. CONSISTENT RESULTS.

**Strong, rigid fixation** superior to bio-absorbable screws

**Radiolucent** for exceptional imaging clarity over metallic screws

**Bio-Inert** properties eliminate inflammatory responses common to absorbable polymers


**No tapping required** facilitated by strong, uniquely angled tip that saves procedural time

**Ideal modulus** of PEEK™ minimizes stress shielding and delivers enhanced bone healing


**Easily revised** through removal or over drilling

## Ordering Information

iFix™ Interference Screws & Disposable Kit	
Model No.	Description
CM-1720	iFix™ Interference Screw - 7mm x 20mm
CM-1725	iFix™ Interference Screw - 7mm x 25mm
CM-1820	iFix™ Interference Screw - 8mm x 20mm
CM-1825	iFix™ Interference Screw - 8mm x 25mm
CM-1830	iFix™ Interference Screw - 8mm x 30mm
CM-1920	iFix™ Interference Screw - 9mm x 20mm
CM-1925	iFix™ Interference Screw - 9mm x 25mm
CM-1930	iFix™ Interference Screw - 9mm x 30mm
CM-1025	iFix™ Interference Screw - 10mm x 25mm
CM-1030	iFix™ Interference Screw - 10mm x 30mm
CM-1130	iFix™ Interference Screw - 11mm x 30mm
CM-1230	iFix™ Interference Screw - 12mm x 30mm
CM-1500	iFix™ Interference Screw Disposable Procedure Kit



7-8mm iFix™ Driver



9-12mm iFix™ Driver



Caution: Federal law (USA) restricts this device to sale only by or on the order of a physician. © Copyright 2010 Cayenne Medical, Inc. All rights reserved. Cayenne Medical is a trademark of Cayenne Medical, Inc. These products are covered by one or more U.S. patents pending. Printed in the USA. iFix™ is a trademark of Cayenne Medical, Inc. PEEK™ is a trademark of Vitrex, PLC.

20322 - Rev. B 2/10

16597 N. 92nd St., Suite 101, Scottsdale, AZ 85260  
cayennemedical.com | P: 888.229.3661

