

The AperFix[®] II System

A COMPLETE
ANATOMIC
SOLUTION



Fast.

Easy, reproducible technique
with rapid deployment

Strong.


Unsurpassed strength utilizing
multidirectional fixation and active
tendon compression at the aperture

Complete.



One system for femoral and tibial
fixation, easy delivery Transtibial or
through the AM Portal and restoration
of native knee biomechanics

PROVEN SYSTEM

AperFix® II Restores Rotational Stability*



AperFix® II Improves Anterior Translation*



* Data upon request.


ENHANCED DESIGN

1479N*



- NEW Innovative solutions for multi-ligament reconstruction techniques
- Unsurpassed graft stiffness replicating the native ACL
- 32% less anterior translation over cortical buttons

1165N*



- NEW cannulated screw & sheath system for precise anatomic positioning
- Circumferential graft compression promoting bone to tendon healing
- Bio-inert, radiolucent PEEK™ polymer providing optimal healing environment

COMPLETE SOLUTIONS

	AperFix® II	Interference Screws	Two Tunnel Double Bundle	Cortical Buttons	Cross-Pin
Aperture Fixation <ul style="list-style-type: none"> ■ Stiffer Graft Construct ■ Decreased chance of tunnel widening 	✓	✓			
Distinct Separation and Positioning of AM & PL Bundles <ul style="list-style-type: none"> ■ Recreates native knee kinematics 	✓		✓		
Fast, Easy Surgical Technique <ul style="list-style-type: none"> ■ Reproducible ■ Saves surgical time 	✓			✓	
Active Circumferential Graft Compression <ul style="list-style-type: none"> ■ Optimal tendon to bone healing environment ■ Decreased chance for windshield wiper effect 	✓				
Easy Removal <ul style="list-style-type: none"> ■ Saves surgical time ■ Maintains original tunnel diameter 	✓				

Ordering Information

AperFix® II Femoral Implant with Inserter


Model No.	Description
CM-2909	9mm x 29mm Femoral Implant with Inserter
CM-2910	10mm x 29mm Femoral Implant with Inserter
CM-2911	11mm x 29mm Femoral Implant with Inserter

AperFix® II Tibial Implant with Driver


Model No.	Description
CM-3008	8mm x 30mm Tibial Implant with Driver
CM-3009C	9mm x 30mm Cannulated Tibial Implant with Driver
CM-3010C	10mm x 30mm Cannulated Tibial Implant with Driver
CM-3011C	11mm x 30mm Cannulated Tibial Implant with Driver

AperFix® Disposable Instruments

Model No.	Description
CM-7014	Calibrated Drill Tipped Guide Wire, 14" x 2.4mm
CM-1500	ACL Disposable Procedure kit



AperFix® II Femoral Implant with Inserter



AperFix® II Tibial Implant with Driver

Caution: Federal law (USA) restricts this device to sale only by or on the order of a physician. Cayenne Medical is a trademark of Cayenne Medical, Inc. These products are covered by one or more U.S. patents pending. All rights reserved. Printed in the USA.

CE

AperFix is a registered trademark of Cayenne Medical, Inc. EndoButton is a trademark of Smith & Nephew. PEEK is a trademark of Victrex, PLC.

11243 - Rev. B