

The iFix[®] System

PEEK-OPTIMA[®] Interference Screw for Bone-Tendon-Bone Fixation

Fast • Strong • Complete



	iFix [®]	Metal Screws	Bio Screws
Strong, Rigid Fixation	✓	✓	
No Tapping Required	✓	✓	
Radiolucent	✓		✓
Bio-Inert	✓	✓	
Ideal Modulus	✓		
Easily Revised	✓		✓



iFix[®] BTB Interference Screw, 10mm x25mm



Cayenne Medical Inc.
16597 North 92nd St, Suite 101 | Scottsdale, AZ 85260
cayennemedical.com | 1-888-229-3661

OPTIMAL PERFORMANCE

CONSISTENT RESULTS

Strong, rigid fixation superior to biocomposite screws.




No tapping required facilitated by strong, uniquely angled tip that saves procedural time.

Radiolucent for exceptional imaging clarity compared to metallic screws.

Bio-Inert properties eliminate inflammatory responses common to absorbable polymers.

Ideal modulus of PEEK-OPTIMA® minimizes stress shielding and delivers enhanced bone healing.

Easily revised through removal or over - drilling.



¹ 20164-C Interference Screw Design Verification & Validation Report

The iFix® BTB Interference Screw is available in 12 sizes. Contact Cayenne Medical for order information.

These products are covered by one or more U.S. patents.



Cayenne Medical Inc.
16597 North 92nd St, Suite 101 | Scottsdale, AZ 85260
cayennemedical.com | 1-888-229-3661